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Type SATURN[®] 50/100/150

for

cryogenic liquid nitrogen and
cryogenic liquid argon

Inhalt	<u>Seite</u>
1 Manufacturer	1
2 Declaration of Conformity	2
4 Description	4
5 Symbols	5
6 Transport/consignment storage/disposal	6
7 Safety	7
8 Vessel drawing	9
9 Technical data	10
10 Preparing the cryovessel for use	12
11 Operation	13
11.3 Transport	15
11.5 Decommissioning	17
11.6 Maintenance	18
11.7 Recurrent/periodic inspections	18
12 Spare parts	19

1 Manufacturer

**Cryotherm GmbH & Co. KG
Euteneuen 4
D - 57548 Kirchen (Sieg)**

**Tel.: + 49- (0)2741/9585-0
Fax: +49- (0)2741/6900
E-mail: info@cryotherm.de
Internet: www.cryotherm.de**

2 Declaration of Conformity

Cryotherm GmbH & Co. KG declares that the conformity assessment and production of the portable cryovessel SATURN[®] 50/100/150 as well as its supply and availability on the market comply with Council Directive 2010/35/EU.

The portable cryovessel is marked with the

0035 marking

3 Marking

on the vessel's identification plate:

- Manufacturer's name and address / type of vessel
- Year of production / year of last recurrent/periodic inspection
- Serial number / production and test guidelines
- Contents
- Maximum operating pressure MOP
- Net weights
- Lowest operating temperature
- Filling material
- Marking in accordance with Directive 2010/35/EU, ADR

On the vessel:

- Label no.11, vessel must be stored and transported in an upright position 

- Hazard warning label no.2, non-inflammable nontoxic gases 

- Label, cryogenic liquid nitrogen
- Label, cryogenic liquid argon
- Label, brief operating instructions

4 Description

Description / intended purpose / intended use

SATURN ® 50/100/150 is a portable cryovessel for the storage and transport / filling and discharge of:

**Cryogenic liquid nitrogen UN no.:1977 or
cryogenic liquid argon UN no.: 1951
under operating pressure.**

The contents must be clearly indicated on the vessel.

The cryovessel must only be operated under the operating conditions named in this operating manual. In particular, the maximum operating pressure indicated must be observed.

The cryovessel is not suitable for outdoor use /operation.

The cryovessel is not suitable for use / operation in ambient temperatures below - 25 °Celsius or above + 40 °Celsius.

Cryotherm GmbH & Co. KG assumes no liability, if the cryovessel is modified or adapted without prior approval of the manufacturer.

5 Symbols

DANGER indicates a high-risk hazard which will result in death or severe injury if not avoided.

WARNING indicates a medium-risk hazard which could result in death or severe injury if not avoided.

CAUTION indicates a low-risk hazard which could result in minor or moderate injury if not avoided.

**6 Transport/consignment
storage/disposal**

- This operating manual is an integral part of the cryovessel and must be kept in a safe place during the operational life span of the cryovessel and given to subsequent operators of the cryovessel.
- After receipt of the cryovessel, the components supplied must be checked for completeness and possible damage.
- If any transport damage is found, immediately inform the transport insurance, the carrier and Cryotherm GmbH & Co. KG.
- In the event of an interim storage of the cryovessel, the vessel must be stored in a dry place at room temperature where it is protected from impact and the penetration of moisture.
- The operator is responsible for the safe disposal of the packaging material.

7 Safety

Cryotherm GmbH & Co. KG recommends that the operator requests the current EC material safety data sheet in accordance with TRGS 220 from the supplier of the substance(s).

cryogenic liquid nitrogen, UN-No.: 1977

cryogenic liquid argon UN no.: 1951 class 2, classification code 3 A

Note:

Any work on these cryovessels must only be carried out by specialist personnel who are familiar with the associated risks and appropriate safety measures.

Possible hazards:

cryogenic liquid nitrogen

cryogenic liquid argon

Cryogenic liquefied gas. When the fluid or large amounts of gas are discharged, a cold fog develops that spreads on the ground. The liquid evaporates quickly as it expands, generating large amounts of gas that can lead to asphyxiation due to the displacement of air. Cold liquid causes frostbite.

Accidental release measures:

In the event of a liquid / gas leakage immediately leave the room, warn others and ensure there is sufficient ventilation. As long as it hasn't been ascertained that the atmosphere is non-hazardous, only enter the area with a self-contained breathing apparatus. In the open, stay on the side facing the wind. Seal off the area.

First aid measures:

Protect yourself and rescue the injured person moving them into fresh air and carefully removing contaminated clothes. Lay down the affected person, keeping them still and warm. If the person is likely to become unconscious put them into the recovery position. In the event of respiratory arrest give mouth-to-mouth or artificial respiration. In the event of respiratory and circulatory arrest start cardiopulmonary resuscitation. Call an emergency physician to the scene of the accident.

In the event of frostbites, hold the person, warming them up with your body heat – do not rub! Do not open blisters, cover wounds with sterile material. If the person's eyes have suffered frost injuries, rinse them carefully opening the eyelid fissure (spasm of the eyelid!)

A legally binding procedure cannot be derived from this information.

The following precautions must be observed:

- The provisions of ADR 2011 regarding transport, transport securing devices etc. must be complied with.
- When setting up in rooms, ensure good ventilation (TRB 610)
- Operation is only permissible to properly trained persons (Technical Rules for Pressure Vessels TRB 700)
- Accident Prevention Regulation - Rule BGR 500 Kap. 2.33 of the German Institution for Statutory Accident Insurance and Prevention (Berufsgenossenschaft).
- TRBS 1111 hazard assessment and safety-related assessment
- TRBS 1201 part 2 – tests in the event of danger from steam and pressure.
- TRBS 2141 danger from steam and pressure – general requirements
- TRG 402 filling systems / operation of filling systems
- TRB 280 General requirements for compressed-gas vessels / operation of compressed-gas vessels
- BGV B6 gases
- Ordinance on Industrial Safety and Health


The regulations and rules for safety and health protection at work issued by the German Institution for Statutory Accident Insurance and Prevention, as amended from time to time, as well as the respective implementation instructions must be observed.

A legally binding procedure cannot be derived from this information.

For safe operation:

- Only operate filling and discharge devices in compliance with the cryovessels' operating conditions.
- Check sealing and functioning of the fittings on a regular basis.
- Use only original spare parts.
- Employ suitable tools.
- Do not operate valves abruptly or jerkily.
Protect lockable spaces from exceeding the maximum operating overpressure by means of a safety valve.
- Adjustment, maintenance and repair work must only be carried out by authorized specialist personnel.
- Do not carry out any mechanical or thermal work on the cryovessel (loss of vacuum).
- Protect safety valves from splash water / lyes.
- Wear personal protective equipment.

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

1. **Protection for fittings**
2. **Pressure build-up valve**
3. **Screw fittings**
(double nipple R 3/8"-3/4-1UNF) for flexible filling hose
4. **Exhaust gas and overflow valve**
5. **Pressure indicator**
6. **Safety valve**
7. **Level indicator**
8. **Combined vacuum seal and safety device**
9. **Transfer siphon with small flange connections (EK siphon)**
10. **Centering and clamping ring**
11. **Pressure rise-control valve *optional**
12. **Pressure drop valve (gas control valve) *optional**

When the display in level indicator is fluctuating,
please adjust the regulating valve in + line so
that the fluctuating calms down

**A complete closing of the regulating valve
leads to incorrect display values**

10 Preparing the cryovessel for use

Any work on this cryovessel must only be carried out by specialist personnel who are familiar with the associated risks and appropriate safety measures.

- Check the components of the cryovessel for completeness and possible damage.
- Check availability of the necessary documents and information.
- Observe possibly necessary permits to work or restrictions.
- Clearly indicate the contents on the vessel with a label.

10.1 Installation

Any work on these cryovessels must only be carried out by specialist personnel who are familiar with the associated risks and appropriate safety measures.

- Control whether all detachable connections are tight.
- Check completeness of the cryovessel, in particular whether the safety devices attached have the correct pressure stage (safety valves) and marking (red mark on the manometer).
- Keep the cryovessel away from electrical devices, sparks, heat sources and open flames.
- Ensure very good ventilation, particularly at the ceiling
- Check tight and reliable fit of the transfer siphon with small flange connection and the holding clamp / clamping ring

11 Operation

Any work on these cryovessels must only be carried out by specialist personnel who are familiar with the associated risks and appropriate safety measures.

Before putting the vessel into service (commissioning) the specialised personnel must be trained regarding the possible dangers, the safety devices and safety measures.

We recommend preparing operating instructions in writing for all activities related to the operation of the cryovessel and making those instructions available to the specialised personnel.


Before each use check the cryovessel for:

- Tight and reliable fit of the transfer siphon with small flange connection and the holding clamp / clamping ring
- Free pass of the safety valve
- Completeness of the fittings and their proper condition

11.1 Filling the vessel

Caution!

- Observe safety notes
- Wear gloves and protective goggles
- Prevent the vessel from rolling away or toppling over and damage


- Close pressure build-up valve 2
- Attach the filling pipeline to filling and discharge valve 4
- Open exhaust gas and overflow valve 5
- Open filling and discharge valve 4
- Finish filling, when liquid leaks from exhaust gas and overflow valve 5
- Close filling and discharge valve 4
- Once the filling and discharge valve is closed, the excess pressure building up inside the vessel due to evaporation is discharged (via overflow valve 13 *optional).

11.2 Discharging the vessel

Caution!

- Observe safety notes
- Wear gloves and protective goggles
- Prevent the vessel from rolling away, toppling over and damage


- Connect withdrawal device to filling and discharge valve 4
- Open pressure build-up valve 2 until the set working pressure (at pressure rise valve 12 * optional) is reached. After that open filling and discharge valve 4.
- After discharging close filling and discharge valve 4.
- Close pressure build-up valve 2
- Open exhaust gas and overflow valve 5 for a possible pressure release of the cryovessel.

11.3 Transport

Cryogenic liquid nitrogen, UN-No.: 1977
Cryogenic liquid argon, UN-No.: 1951

- Hazard warning label 2.2, non-inflammable non-toxic gases
- Hazard warning label 11, two black arrows
- ADR / GGVSEB / GGVSee / class 2, classification code 3 A
- ADR hazard number 22
- Packaging instructions P 203
- Ensure the vehicle driver is aware of the possible hazards of the load and knows what to do in the event of an accident or emergency.
- Avoid transport on vehicles where the loading space is not separated from the driver's compartment.
- Only transport full and empty vessel with closed and tight valves.
- Ensure that the cryovessels are standing upright and prevented from shifting or falling over.

- During the transport of the cryovessel pressure build-up valve 2, filling and discharge valve 4 and exhaust gas and overflow valve 5 must be closed.
- We recommend transporting the cryovessel when full and with a vessel pressure of max. 1 bar.

11.4 Malfunction

Immediately decommission vessel if

- The holding clamp / the clamping ring at the transfer siphon with small flange connection is damaged or unusable
- Fittings are leaking
- The safety valve is blowing off intensively
- The outer cylinder shows condensation or icing, this indicates vacuum loss

Vessels with a vacuum loss are unusable.

Return to manufacturer for examination / repair.

Should you have any queries, please indicate

- Type of vessel
- Fabrication number
- Year of construction

In the event of operating troubles always inform the manufacturer.

Manufacturer **Cryotherm GmbH & Co. KG**
Euteneuen 4
57549 Kirchen (Sieg)
Tel.: 02741 / 9585-0
Fax: 02741 / 6900

11.5 Decommissioning

- If the vessel is taken out of service, it must be emptied completely, warmed up and stored with a slight overpressure of gas in order to prevent the condensation of moisture.
- Before dismantling the vessel, check whether it is in an unpressurised condition.
- If the cryovessel is permanently decommissioned, the components must be recycled / disposed
- of separately

11.6 Maintenance

- **Maintenance work on the cryovessel must only be carried out by trained specialist personnel with the appropriate equipment/tools.**
- Maintenance work on the cryovessel should include a check of
 - The proper condition of the holding clamp / the clamping ring on the transfer siphon with small flange connection
 - Icing of connections, couplings or screw connections
 - Complete labelling
 - Tightness of connections, couplings or screw connections
 - Functioning, opening pressure and sealing of the safety valve
 - Functioning of the manometer
 - Control of the swivel castors
- Cryotherm recommends that maintenance work on the cryovessel is carried out every 2 years.

11.7 Recurrent/periodic inspections

- **Test intervals acc. to EN 1251-3 operating requirement every 10 years**
- **In addition, according to EN 1251-3 the response of all safety valves must be checked by a competent authority or body at intervals of no longer than 5 years.**

12 Spare parts

Item	Name	Item number
	SATURN®	
7	safety valve Herose 3.0 bar	78214249
6	manometer 0 to 4 bar WIKA	0640556
2	ball valve 3/8"	0346570
8	level indicator WIKA	78202720
13	overflow valve, opening Herose 1 – 4 bar *optional	78213670
4	positive displacement valve 3/8" Rego	78213908
5	positive displacement valve 1/2" Rego	78213459
12	pressure build-up regulator * option	0366006
	filling hose	79229957
	phase separator	79407217
11	outer centering ring DN 50	0321303
11.1	KF clamping ring DN 50	0792277
	ADR – marking	
	nitrogen label	78400571
	argon label	78202942
	ADR label no. 2	0358193
	ADR label ↑↑ no. 11	0356199
	operating manual	78208074
	label brief operating instructions	79406614

Cryotherm

Cryotherm GmbH & Co. KG certified according to DIN EN ISO 9001:2008
Article - No. :78208074 - 1019
Subject to change
© Cryotherm GmbH & Co. KG
® registered trademark

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Cryotherm GmbH & Co. KG Deutschland